

## Step 1:



Begin by flattening the Flo-Form fluidized positioner.

## Step 2:



Position baby near top end of the positioner and mold Flo-Form toward the baby.

## Step 3:



With "Mothers Hand" on the baby, use opposite hand to mold the far side of the positioner toward the baby.

## Step 4:



Repeat Step 3, molding the near side of the positioner toward the baby.

## Step 5:



Mold top of positioner to head against hand.

## Step 6:



Mold bottom of positioner to feet against hand and place tether straps where appropriate.



## For Proning:

Build prone roll from hips to top of head by pushing down with fingers along side of body. Compress positioner below umbilicus and knees. A small utility positioner can be used directly on top of full body positioner as a prone roll.



## For Side-Lying:

Mold Flo-Form against the outer perimeter of the baby's body, being sure to avoid any lines and tubing.

Note: Flo-Form neonatal positioners are recommended for single patient use.

v.1.1 5/25/22