

Lateral Turn Hold



Hold positioner horizontally by the corners so the beads fall to the bottom



Roll patient onto their side



Tuck the leading edge of the positioner under the patient's back



Roll patient back onto the positioner

Floating the Heels



Place flattened positioner under the calf/calves to offload the heels for 1-2" of lift



Stand positioner up on its long edge and place under the calf/calves for 3-4" of lift

Offloading the Sacrum



Place flattened positioner under the patient's thighs to partially offload the sacrum

DO:

- Clean with hospital-approved disinfectant
- Use multiple positioners for larger patients
- Position with filter facing any direction

DON'T:

- Stuff or fold the positioner under the patient
- Use a product if it is punctured or damaged