

## What is a fluidized positioner?

Flo-Form fluidized technology allows the positioners to be molded and contoured to each infant's changing size, medical condition, and equipment. This individualized positioning provides comfort and support to vulnerable infants.

## What are the thermal regulating properties of fluidized positioners?

Fluidized positioners have a very low thermal mass and will assume the ambient environmental temperature very quickly. For babies in a heated environment, pre-warm positioner under the radiant warmer or in the incubator before use.

## What is the difference between the grey and blue positioners?

Functionally they allow for the same positioning and micro-contouring, but the blue material is slightly softer to the touch and has two layers of material for additional anti-shear benefits.

## How do you clean a fluidized positioner?

Flatten the positioner and smooth wrinkles from the surface. Wipe all sides of the positioner with a facility-approved disinfectant. Allow it to dry before placing in a clean cover or pillowcase. Cloth covers can be laundered in very hot water and dried on low heat for reuse.

## Why use fluidized positioners?

Fluidized positioners encourage consistent flexed / midline containment, allowing for optimal neuromotor development and relaxation. This also has a positive impact on physiologic stability, skin integrity, and sleep.

## How does a fluidized positioner compare to other positioning aids?

The versatility and effectiveness of fluidized positioners can help reduce the need for multiple positioning products (gels, blankets, foam, and cloth positioning aids). A single fluidized positioner can typically address a neonates complete positioning needs.

## Can you X-ray through a fluidized positioner?

Spry does not recommend that you x-ray through a fluidized positioner.

## When should you dispose of a fluidized positioner?

A positioner should be disposed of once it is not able to maintain its position when molded.

## What happens if the positioner is torn or punctured?

Simply replace the positioner. Do not attempt to repair the product.